



# GABE FAUBERT

Midwestern US & Eastern Canada Regional Manager  
R-Biopharm Inc.

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission



# EVALUATING PROBLEMATIC MATRICES

Gabe Faubert

Eastern Canada Regional Manager, R-Biopharm, Inc.

9-27-16



THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada




Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# R-BIOPHARM



THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



# R-BIOPHARM & GLUTEN, A BRIEF HISTORY

<b>1988</b>	Foundation of Rohm GmbH (later changed to R-Biopharm)
<b>1998</b>	Foundation of R-Biopharm, Inc (USA & Canada)
<b>1999</b>	R5 Sandwich ELISA prototype created
<b>2001</b>	Launch of R-Biopharm RIDASCREEN® Gliadin R5 ELISA test kit
<b>2002</b>	Launch of R-Biopharm RIDA® Quick Gliadin R5 lateral flow test kit
<b>2006</b>	RIDASCREEN® Gliadin R5 ELISA received AOAC-RI approval
<b>2007</b>	Launch of R-Biopharm R5 Gliadin Competitive test kit (Gen1)
<b>2008</b>	Codex Alimentarius: 20ppm threshold and R5 sandwich as Type 1 Method
<b>2009</b>	Launch of R-Biopharm Surefood PCR Gluten
<b>2011</b>	Launch of updated R-Biopharm R5 Gliadin Competitive test kit (Gen2)
<b>2012</b>	RIDASCREEN® Gliadin R5 ELISA is adopted as an AOAC Official Method (2012.01)
<b>2015</b>	RIDASCREEN® Gliadin Competitive is adopted as an AOAC Official Method (2015.05)
<b>2016</b>	RIDA® QUICK Gliadin R5 Lateral Flow is adopted as an AOAC Official Method (2015.16)

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# CASE 1 - OATS

**Question: Does the R5 antibody react with Oats?**

**Origin: Oat samples occasionally testing positive**

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# CASE 1 - OATS

## Experiment

*Gather a variety of certified gluten free oat samples*

*Analyze them with multiple antibodies*

*Confirm positive results with PCR*

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# CASE 1 - OATS

Oat variety	Antibody				PCR
	R5 Sandwich ELISA (R7001)	Antibody 2 Sandwich ELISA	Antibody 3 Sandwich ELISA	Antibody 4 Sandwich ELISA	PCR (on wheat, rye and barley - multiplex)
Ivory	positive	positive	positive	positive	positive
Scorpio	positive	positive	positive	positive	positive
Zorro	negative	negative	negative	negative	Not tested
Poseidon	negative	negative	negative	negative	Not tested
Symphony	positive	positive	positive	positive	positive
Oberon	negative	negative	negative	negative	Not tested
Aragon	negative	negative	negative	negative	Not tested
Husky	negative	negative	negative	negative	Not tested
Flämingslord	negative	negative	negative	negative	Not tested
Lutz	negative	negative	negative	negative	Not tested

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



# CASE 1 - OATS

## Conclusion

Pure gluten free oats do not generate a positive response with the R5 antibody

## *Notes:*

The 3 positive varieties were reacquired at a later date and tested negative

This work was done in 2013 and more varieties have since been tested and verified

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 



# CASE 2 - HORSERADISH

**Question: Does horseradish cause false positives due to horseradish peroxidase (HRP)?**

**Origin: Certain horseradish powders were tested positive by ELISA**

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# CASE 2 - HORSERADISH

## Experiment:

*Test “positive” horseradish powder with Antibody and PCR methods to determine contamination...or false positive due to HRP*

## **Samples to be analyzed:**

- #1 ‘negative’ horseradish powder (*customer supplied*)
- #2 ‘positive’ horseradish powder (*customer supplied*)
- #3 ‘very positive’ horseradish powder (*customer supplied*)
- #4 Raw horseradish (*purchased from grocery store*)

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# CASE 2 - HORSERADISH

Sample	Method			
	RIDASCREEN®Gliadin (Ethanol prep)	RIDASCREEN®Gliadin (Cocktail prep)	RIDASCREEN® QUICK Gliadin (ethanol prep)	SureFood® ALLERGEN ID Gluten PCR
#1 negative HR powder	<5ppm gluten	<5ppm gluten	<5ppm gluten	Negative
#2 positive HR powder	75ppm	>80ppm	~80ppm (serial dilution)	Positive
#3 very positive HR powder	~1500ppm	>80ppm	>160ppm (serial dilution)	Positive
#4 raw HR	<5ppm gluten	<5ppm gluten	<5ppm gluten	Negative

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



# CASE 2 - HORSERADISH

## Conclusion

Horseradish peroxidase contained in horseradish does not cause false positives in either RIDASCREEN® Gliadin ELISA or RIDASCREEN® QUICK Gliadin Lateral Flow

## **Notes:**

For sample #2 it is possible that accidental contamination was to blame, however adulteration was likely the cause of the contamination in sample #3

AFTER the science proved contamination, the ingredient broker investigated his supplier in Asia...

And our customer switched brokers

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



# CASE 3 - PROTEASE

**Question: Does a protease concentrate cause false negatives in sandwich ELISA and false positives in competitive ELISA?**

**Origin: An enzyme manufacturer wanted to ensure no carryover from a potential contamination sources during the manufacturing process. Their analysis indicated interference when testing the enzyme concentrate.**

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# CASE 3 - PROTEASE

## Experiment:

*Conduct a series of spike & recovery analyses to determine root cause*

Due to the known interference of enzyme concentrates it is necessary to include a denaturing step in the sample prep.

All experiments were tested with a sandwich ELISA (RIDASCREEN®FAST Gliadin) and competitive ELISA (RIDASCREEN®Gliadin competitive)

*Experiment #1      Test enzyme with normal sample prep*

*Experiment #2      Test enzyme with denaturing sample prep*

*Experiment #3      Spike enzyme before denaturing sample prep*

*Experiment #4      Spike enzyme after denaturing sample prep*

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



# CASE 3 - PROTEASE

Experiment	Method			
	<u>RIDASCREEN® FAST Gliadin</u>		<u>RIDASCREEN® Gliadin competitive</u>	
	Result	% recovery	Result	% recovery
Normal prep	<10 ppm		132.4 ppm	
Denatured prep	<10 ppm		43.1 ppm	
Spike <i>before</i> denature	<10 ppm		195.6 ppm	316%
Spike <i>after</i> denature	23.5 ppm	118%	75.7 ppm	163%
Spike control	35.2	136%	37.5 ppm	187%

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



# CASE 3 - PROTEASE

## Conclusion

The untreated enzymes will cause false negatives in a sandwich ELISA and false positives in a competitive ELISA

It is theorized that the heating process accelerated the enzyme activity causing it to alter the spike before being denatured

## **Note:**

Results from RIDASCREEN® Gliadin competitive indicated that there were still allergenic peptide fragments present in the sample, thus confirming contamination

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:





# CONCLUSION

Through 15+ years of servicing the gluten free industry, R-Biopharm has had the opportunity to evaluate many unique samples.

This experience gives us great confidence that we will find a solution for any potentially problematic sample matrix.

THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:



Agriculture and  
Agri-Food Canada



Canadian Grain  
Commission

Growing Forward 2 



Pulse Canada 

# Thank You



THE 2016 GLUTEN-FREE STAKEHOLDER UPDATE & PLANNING SESSION PRESENTED IN COOPERATION WITH:

